

SKID-E Evacuation Chair



The most important innovation is the use of passive slides which glide over the staircase with a high standard of safety and without requiring any effort. An attentive study of the chair underlines a high esthetical value with a new concept of style.

The rugged look is a mix of the technical details found in the unique form of the frame or the clearly original position of the seat. Precious aesthetic details are offered by the angled parts, the inclination of the telescopic backrest, the headrest adjustable in height and the leaning transport handle which is distinctively the continuation of the outline.

Design, innovation and maximum safety in use characterize the new Skid- Series E. An absolute novelty in the world of evacuation and transport chairs.

A peculiarity of the Skid-E is that there is no support to maintain it in the vertical position. An attentive study of new geometrical solutions offers the device complete stability even in this mode.

Skid-E is the product of a long term project which has made it extremely light and versatile and asks for no compromise regarding the conditions in which it can be used. The project has also studied solutions permitting maximum manoeuvrability in confined spaces and on staircases making staircase transport and building exit a simple operation. The innovative sliding system allows the operator to control descent of the staircases without having to support any weight.

Manufactured in aluminium with a special polished finish. Seat covers are in Spentex®.

Front pivoting wheels with break Ø 100 x 32

Fixed rear wheels Ø 200 x 50

Skid-E is supplied with two polypropylene belts 50 mm.

Technical details

| | |
|----------------|-----------|
| Material | Aluminium |
| Open height | 1600 mm |
| Folded height | 1140 mm |
| Width | 530 mm |
| Width folded | 330 mm |
| length | 900 mm |
| Top finish | Polished |
| Weight | 12,7Kg |
| Loading weight | 150Kg |

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
E-mail: info@spencer.it Web: http://www.spencer.it

This document has been automatically generated by www.spencer.it on 2009/09/11 09:09:10, and is referred to the last production lot. The informations furnished have been last updated on . To verify the correctness of this document, please contact our Quality Dept. at quality@spencer.it.

SK20001E

Skid-E Silver

