

**SKID-E READY**  
Evacuation Chair



Skid-E Ready is characterised by the new fork and tubular frame, which enhances the technological design of the assembly and transfer of the castor wheels at the front of the chair, which is divided into two elements of very large surface tilted back and protected by two lateral conveyors that integrate the footrest.

Made of aluminum, the Skid-E Ready, comes with the upper handle fixed to be always ready for use.

The chair Skid-E Ready comes with two 50 mm polypropylene belts.

**TECHNICAL DATA**

Height	1580 mm
Width	540 mm
Length	930 mm
Thickness folded	320 mm
Weight	12,5 kg
Loading capacity	150 kg
Surface treatment	Painting
Anterior wheels	n°2 Ø 100 mm castor wheels with brake
Posteriori wheels	n°2 Ø 200mm
Canvas	Spentex ®
Canvas color	Black
Frame	Alluminum
Frame color	Yellow
Belts	50mm polypropylene
Latex free	

SK20101C

SKID-E READY



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>